



medical
magnesium

mm.CS

Compression Screw


Bioabsorption. Stability of metal.

As a metallic and resorbable double-threaded screw, the **mm.CS** fixes bone fragments securely and primarily stable in order to ensure a controlled healing process.

The **mm.CS** is available in three diameters (2.8, 3.5 and 5.0 mm) and various lengths to ensure correct sizing for all intended indications. Please see following page for exact configurations and corresponding order information.

Intended purpose

mm.CS is a bioabsorbable compression screw used for reconstructive procedures to reinstate bone continuity during osteosynthesis of fractured bones and osteotomies of small bones, to refixate smaller cortical, cancellous and osteochondral fragments or avulsion fractures, as well as arthrodesis of smaller joints during therapeutic open surgical, arthroscopic and minimally invasive orthopedic and/or trauma surgery.





Size 2.8

Length	Distal Thread length	Article number
10	4	011.281004
12	4	011.281204
14	4	011.281404
16	4	011.281604
18	5	011.281805
20	5	011.282005
20	8	011.282008
22	6	011.282206
22	10	011.282210
24	6	011.282406
24	10	011.282410
26	7	011.282607
26	11	011.282611
28	7	011.282807
28	11	011.282811
30	8	011.283008
30	12	011.283012
32	8	011.283208
32	12	011.283212
34	9	011.283409
34	14	011.283414

Size 3.5


Length	Distal Thread length	Article number
10	4	011.351004
12	4	011.351204
14	4	011.351404
16	4	011.351604
18	5	011.351805
20	5	011.352005
22	6	011.352206
24	6	011.352406
24	10	011.352410
26	7	011.352607
26	11	011.352611
28	7	011.352807
28	11	011.352811
30	8	011.353008
30	12	011.353012
32	8	011.353208
32	12	011.353212
34	9	011.353409
34	14	011.353414
36	9	011.353609
36	14	011.353614
38	10	011.353810
38	15	011.351815
40	10	011.354010
40	15	011.354015

Size 5.0

Length	Distal Thread length	Article number
12	4	011.501204
14	4	011.501404
16	4	011.501604
16	6	011.501606
18	5	011.501805
20	5	011.502005
22	6	011.502206
24	6	011.502406
24	10	011.502410
26	7	011.502607
26	11	011.502611
28	7	011.502807
28	11	011.502811
30	8	011.503008
30	12	011.503012
32	8	011.503208
32	12	011.503212
34	9	011.503409
34	14	011.503414
36	9	011.503609
36	14	011.503614
38	10	011.503810
38	15	011.503815
40	10	011.504010
40	15	011.504015
42	11	011.504211
42	17	011.504217
44	11	011.504111
44	17	011.504417
46	12	011.504612
46	18	011.504618
48	12	011.504812
48	18	011.504818
50	14	011.505014
50	20	011.505020
55	14	011.505514
55	22	011.505522
60	15	011.506015
60	24	011.506024

Size 2.8

Cannula 1.1
K-Wire 1.0
Drive T8



Size 3.5

Cannula 1.2
K-Wire 1.1
Drive T8



Size 5.0

Cannula 1.7
K-Wire 1.6
Drive T15

