



medical
magnesium

mm.CS

Compression Screw


Bioabsorption. Stability of metal.

As a metallic and resorbable double-threaded screw, the **mm.CS** fixes bone fragments securely and primarily stable in order to ensure a controlled healing process.

The **mm.CS** will be available in three diameters (2.8, 3.5 and 5.0 mm) and various lengths to ensure correct sizing for all intended indications. Please see following page for exact configurations and corresponding order information.

Intended purpose

mm.CS is a bioabsorbable compression screw used for reconstructive procedures to reinstate bone continuity during osteosynthesis of fractured bones and osteotomies of small bones, to refixate smaller cortical, cancellous and osteochondral fragments or avulsion fractures, as well as arthrodesis of smaller joints during therapeutic open surgical, arthroscopic and minimally invasive orthopedic and/or trauma surgery.



Size 2.8

Length	Distal Thread length	Article number
10	4	011.281004
12	4	011.281204
14	4	011.281404
16	4	011.281604
18	5	011.281805
20	5	011.282005
22	6	011.282206
24	6	011.282406
26	7	011.282607
28	7	011.282807
30	8	011.283008
32	8	011.283208
34	9	011.283409

Size 3.5


Length	Distal Thread length	Article number
10	4	011.351004
12	4	011.351204
14	4	011.351404
16	4	011.351604
18	5	011.351805
20	5	011.352005
22	6	011.352206
24	6	011.352406
26	7	011.352607
28	7	011.352807
30	8	011.353008
32	8	011.353208
34	9	011.353409
36	9	011.353609
38	10	011.353810
40	10	011.354010

Size 5.0

Length	Distal Thread length	Article number
32	8	011.503208
34	9	011.503409
36	9	011.503609
38	10	011.503810
40	10	011.504010
42	11	011.504211
44	11	011.504411
46	12	011.504612
48	12	011.504812
50	14	011.505014

Size 2.8

Cannula 1.1
K-Wire 1.0
Drive T8



Size 3.5

Cannula 1.2
K-Wire 1.1
Drive T8



Size 5.0

Cannula 1.7
K-Wire 1.6
Drive T15

